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Références

- California Department of Public Health. Frequently Asked Questions about the Use of Stockpiled N95 Filtering Facepiece Respirators for Protection from COVID-19 Beyond the Manufacturer-Designated Shelf Life (mise à jour 10 mars, 2020) <https://www.cdph.ca.gov/Programs/CID/DCDC/Pages/FAQ-N95.aspx>
- California Department of Public Health. Cal/OSHA Interim Guidance on Coronavirus for Health Care Facilities: Efficient Use of Respirator Supplies. (mars 2020) <https://www.dir.ca.gov/dosh/Use-of-Respirator-Supplies.html>
- Centers for Disease Control and Prevention (CDC). Release of Stockpiled N95 Filtering Facepiece Respirators Beyond the Manufacturer-Designated Shelf Life: Considerations for the COVID-19 Response. (28 février 2020) <https://www.cdc.gov/coronavirus/2019-ncov/release-stockpiled-N95.html#f1>
- Centers for Disease Control and Prevention (CDC). Strategies for optimizing use of N95 masks. (mise à jour 29 février 2020) https://www.cdc.gov/coronavirus/2019-ncov/hcp/respirators-strategy/index.html?CDC_AA_refVal=https%3A%2F%2Fwww.cdc.gov%2Fcoronavirus%2F2019-ncov%2Fhcp%2Frespirator-supply-strategies.html
- Centers for Disease Control and Prevention (CDC). National Personal Protective Assessment Laboratory. Personal Protective Equipment Conformity Assessment Studies and Evaluations (PPE CASE) Reports: Performance of Stockpiled Air-Purifying Respirators. (1er mars 2020) <https://www.cdc.gov/niosh/npptl/ppecase.html>

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- National Institute for Occupational Safety and Health, Centers for Disease Control and Prevention. [Recommended Guidance for Extended Use and Limited Reuse of N95 Filtering Facepiece Respirators in Healthcare Settings](#)
- National Institute for Occupational Safety and Health, Centers for Disease Control and Prevention. [Release of Stockpile N95 Filtering Facepiece Respirators Beyond the Manufacturer-Designated Shelf Life: Consideration for the COVID-19 Response.](#)
- Centers for Disease Control and Prevention. [Strategies for Optimizing the Supply of N95 Respirators](#)
- California Department of Public Health. [Implementing Respiratory Protection Programs in Hospitals: A Guide for Respirator Program Administrators](#)
- California Department of Public Health. [Respirator Use in Health Care – a Toolkit for Program Administrators](#)
- California Department of Public Health. [Respirator Selection Guide for Aerosol Transmissible Diseases](#)
- California Department of Public Health. [Cal/OSHA's Aerosol Transmissible Disease Standards and Local Health Departments](#)
- OSHA. [Hospital Respiratory Protection Program Toolkit: Resources for Respirator Program Administrators](#)
- Cal/OSHA. [The California Workplace Guide to Aerosol Transmissible Diseases](#)
- Cal/OSHA. [Respiratory Protection Factsheet](#)
- Cal/OSHA. [Respiratory Protection in the Workplace: A Practical Guide for Small-Business Employers](#)
- Cal/OSHA. [Interim Guidance for Protecting Health Care Workers from Exposure to 2019 Novel Coronavirus](#)
- Cal/OSHA. Aerosol Transmissible Diseases Standard, [title 8 section 5199](#)

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